



LuminolPenTM Transforms Your Western Blots

LuminolPen™, HRP System is a specialized tool for Western Blot experiments, designed to take notes on blots and to bridge the gap between pre-stained protein markers and chemiluminescent signals.



- Compatible with PVDF/NC
- Works with film/imaging systems
- Improves data management



signal on blots

LuminolPen™, HRP System

#LH03-50 1 pen for 1,000 membranes #LH03-10 1 pen for 100 membranes














Easily differentiate between similar image data

The unique ink of LuminolPen™ is visible in blue and luminescent under chemiluminescent substrates, allowing researchers to annotate protein ladders before detection. This enables simultaneous visualization of markers and protein signals in one image, enhancing accuracy and efficiency in Western blot analysis.





Marker positions

By anotating each blots, it's easy to differentiate between similar image data.

- Pre-stained marker positions
 Sample names
 Experimental dates
- Lane numbers
- Antibodies used Investigator's name

Application of LuminolPen™, HRP System




Protein gel electrophoresis

- Include pre-stained molecular weight markers alongside your samples during electrophoresis.



Protein transer

- After transferring proteins to membrane, gently remove excess liquid from blot with filter paper. Keep the blot slightly damp, avoid over-drying.



Annotate with LuminolPen™

- Use LuminolPen™ to mark visible pre-stained protein bands on your membrane. Add notes or labels on the membrane margins in needed.



Complete 4 Western Blot

- Proceed with blocking and subsequent steps of your Western Blot protocol.



Visualize results

- Use film or imaging system to capture signals of the protien targert and the marks of LuminolPen™.

Reference

- 1. Rahman Md Matiur, et al. Animals 11.9 (2021): 2506. PMID: 34573471.
- 2. Rahman Md Matiur, et al. Membranes 13.2 (2023): 185. PMID: 36837688.
- 3. Nguyen, Tan Phat, et al. Fish & Shellfish Immunology 139 (2023): 108902. PMID: 37330026.



