

Western Blot Transfer Buffer (pH 8.3)

WTB10-1L/WTB500/WTB500P

V2.1

For Research Use Only

Introduction

Western Blot Transfer Buffer (pH 8.3) is used for general SDS-PAGE gels running. When diluted/reconstituted to a 1X working solution, it contains 25mM Tris, 192mM Glycine.

Product Components

Western Blot Transfer Buffer (10X) (WTB10-1L)

Western Blot Transfer Buffer (10X)	500 mL	2 bottles
User's manual		

Western Blot Transfer Buffer, powder (WTB500P)

Western Blot Transfer Buffer	WTB500	500 mL dilute	20 packs
User's manual			

Safety Information

Please wear gloves, lab coat and goggles while operating. Prevent contact product directly. In case of contacting, wash with large amount of water.

Storage

Store **Western Blot Transfer Buffer (pH 8.3)** at 2-8 °C or room temperature. Expiration date is labeled on the bottle/pack or box.

Reconstitution

Reconstitute one part of **WTB10-1L** with nine parts of ddH₂O or one pack of **WTB500** with 500 mL of ddH₂O to make 1X Western Blot Transfer Buffer. If methanol is necessary, please adjust **methanol about 10-20% (v/v)**. The pH of the buffer should be 8.3 and no pH adjustment is required.

Related Visual Protein Products

SDS-PAGE Running Buffer (10X)	RB10-1L	1 L
SDS-PAGE Running Buffer, powder	RB500P	1 set